

Get the surgical table that best suits your needs

Take advantage of ISI's free evaluation program

Free Evaluations ISI offers a choice of tops, legs and configurations

We want you to get the surgical table that best suits your needs. That's why we have a program for qualified facilities to evaluate ISI surgical tables for a 30 day, no-charge evaluation in the O.R. Begin your free evaluation by calling ISI at 1-800-537-1945 or by visiting www.isisurgery.com/ISIOR/SpecialtyTables.htm.

Design a custom table by specifying the features you want

ISI offers the option of creating your own table. Table top material and size, type of leg, mounting bracket and pad – you specify what you want. Note: Custom tables are not available for evaluation.

Standard Feature	Advanced Universal K	Universal K	Carter Hand Surgery Table	Self-Supporting Arm Table	Ultra-Light	Free-Standing Arm Table
Lexan Top		●		●		●
Carbonfiber Top	●				●	
Polycarbonate Top			●			
Weight	16 lbs	20 lbs	18 lbs	18 lbs	9 lbs	25 lbs
Telescoping Leg	●	●	●	optional	optional	●
Side Rail Mounting	●	●	●	●	●	
Fluoroscopy	●	●	●	●	●	●
Traction	optional	optional	standard	optional	optional	optional
Under-Pad Mounting Piece	●	●			●	●
2" & 4" Foam Pads	optional	optional	1", 2", 3" pads optional*	optional	optional	optional
Height Adjustable Support Leg	●	●	●	optional	optional	●
Free Evaluation	●	●	●	●	●	●

*The Carter Hand Surgery Table comes standard with a one-inch foam pad.