

ISI Brush Nest – Brush set and organizer for instrument cleaning brushes

The ISI Brush Nest is the wall organizer with a complete set of 18 instrument cleaning brushes in a variety of sizes, diameters and composition. Made of tough, durable plastic with a convenient perforated shelf, the ISI Brush Nest is the ultimate organizer. Technicians can get to work right away with the right brushes for the task.

The long, thin brushes are perfect for cannulas and tubing. The smallest brushes are the answer for micros. The brass and stainless steel brushes make quick work of box locks.

Individual packages of replacement brushes are available for quick delivery.



A DIVISION OF INSTRUMENT SPECIALISTS, INC.

1-800-537-1945 www.isisurgery.com

| Stock No. | Color Code | Description of Overall Length and Diameter | # Per Pkg | Price Per Pkg |
|-------------------------|------------|---|-----------|---------------|
| BN-9600 | | Brush Nest set – the wall hanger with 18 different brushes | 1 | \$235.00 |
| BN-9650 | | Re-order set of one of each standard brush | 18 | 62.00 |
| Standard Brushes | | | | |
| BN-130-1 | White | 21" tube brush, 20 French, 1/4" diameter with 4" bristle length | 3 | \$18.50/pkg |
| BN-130-3 | Black | 21" tube brush, 30 French, 3/8" diameter with 4" bristle length | 3 | 18.50/pkg |
| BN-56740 | Blue | 20" tube brush, 24 French, 5/16" diameter with 3-1/2" bristle length | 3 | 17.25/pkg |
| BN-660-1 | Red | 6-1/4" tube brush, 9 French, 1/8" diameter with 2-1/4" bristle length | 6 | 16.25/pkg |
| BN-660-5 | Blue | 6-1/4" tube brush, 24 French, 5/16" diameter with 2-1/4" bristle length | 6 | 16.25/pkg |
| BN-661 | Yellow | 5" tube brush, 21 French, 9/32" diameter with 1-1/2" bristle length | 6 | 15.75/pkg |
| BN-662 | Red | 5" tube brush, 9 French, 1/8" diameter with 1-1/2" bristle length | 6 | 15.75/pkg |
| BN-663 | Green | 5" tube brush, 14 French, 3/16" diameter with 1-1/2" bristle length | 6 | 15.75/pkg |
| BN-670 | White | 6-1/4" tube brush, 20 French, 1/4" diameter with 2" bristle length | 6 | 16.25/pkg |
| BN-671 | Black | 6-1/4" tube brush, 30 French, 3/8" diameter with 2" bristle length | 6 | 16.25/pkg |
| BN-940-1 | White | 12" tube brush, 20 French, 1/4" diameter with 4" bristle length | 3 | 13.85/pkg |
| BN-940-3 | Black | 12" tube brush, 30 French, 3/8" diameter with 4" bristle length | 3 | 13.85/pkg |
| BN-SX3198 | Gray | 16" tube brush, 34 French, 1/2" diameter with 3" bristle length | 3 | 15.25/pkg |
| BN-SX3923 | Brown | 13" tube brush, 34 French, 1/2" diameter with 3" bristle length | 3 | 14.85/pkg |
| BN-069 | | Nylon brush – tooth brush style | 12 | 25.00/pkg |
| BN-83112 | | Brass brush – tooth brush style | 12 | 25.00/pkg |
| BN-83111 | | Stainless steel brush – tooth brush style | 12 | 25.00/pkg |
| BN-3000 | | Nylon nail scrub brush | 3 | 9.50/pkg |
| Optional Brushes | | | | |
| BN-4100 | | Pipe cleaners; 6" long, 3/16" diameter | 36 | \$ 3.00/pkg |
| BN-4200 | | Pipe cleaners; 12" long, 3/16" diameter | 24 | 4.00/pkg |
| BN-4300 | | Pipe cleaners; 12" long, 3/8" diameter | 24 | 7.50/pkg |
| B-06213 | | Small trach brushes 3-3/4"; fits tube sizes 000-5, 1-7/8" long, 1/2" diameter | 3 | 12.00/pkg |
| B-06214 | | Large trach brushes 8"; fits tube sizes 6-12; 3-3/4" long, 1/2" diameter | 3 | 14.00/pkg |
| BN-130-18C | Red | 21" tube brush, 9 French, 1/8" diameter with 2" bristle length | 3 | 21.50/pkg |
| BN-56737 | Blue | 14" tube brush, 3/4" diameter with 3-1/2" bristle length | 3 | 15.25/pkg |
| BN-56738 | Yellow | 14" tube brush, 7/8" diameter with 4" bristle length | 3 | 15.25/pkg |
| BN-56739 | Blue | 8" tube brush, 3/4" diameter with 3-1/3" bristle length | 6 | 19.85/pkg |
| BN-660-316 | Green | 6-1/4" tube brush, 14 French, 3/16" diameter with 2-1/4" brush length | 6 | 16.25/pkg |
| BN-660-14 | White | 6-1/4" tube brush, 20 French, 1/4" diameter with 2-1/4" brush length | 6 | 16.25/pkg |

Minimum Order: \$25

Freight: F.O.B. Origin. Freight prepaid and added to invoice.

Shipping: UPS Ground Service, unless otherwise requested.

Terms: Net 30 days

Warranty: The ISI Brush Nest wall organizer is guaranteed for 90 days against manufacturer's defects. The brushes are not covered under the 90 day warranty.


ISI Micro Nest – Brush set and organizer for micro brushes

The ISI Micro Nest is the wall-mounted organizer with a complete set of micro brushes in a variety of sizes and diameters. The ISI Micro Nest makes it easier for technicians to select the right brush for the cleaning task.

ISI Micro Brushes are designed specifically for those hard-to-clean instruments such as

cannulas. ISI Micro Brushes are available in a variety of French diameters. Made with stainless steel handles and nylon bristles, the ISI Micro Brushes will last through countless cleanings.

Brush packages may be ordered separately.



ISI Medical & O.R. Products

A DIVISION OF INSTRUMENT SPECIALISTS, INC.

1-800-537-1945 www.isisurgery.com

| Stock No. | Color Code | Description of Overall Length and Diameter | # Per Pkg | Price Per Pkg |
|----------------------------|------------|--|-----------|---------------|
| M-2004 | n/a | ISI Micro Nest set – the wall organizer with 8 micro brushes (one each listed below) | 1 | \$85.00 |
| M-2800 | n/a | Micro Nest Wall Organizer – the organizer without brushes | 1 | 25.00 |
| M-2650 | n/a | Re-order set of one of each micro brush | 9 | 62.00 |
| Micro Nest Brushes: | | | | |
| B-4FR10 | Gray | Micro Brush, 4 French diameter, 10" overall length | 3 | \$23.50/pkg |
| B-5FR10 | Yellow | Micro Brush, 5 French diameter, 10" overall length | 3 | 23.50/pkg |
| B-6FR10 | Blue | Micro Brush, 6 French diameter, 10 overall length | 3 | 23.00/pkg |
| B-7FR10 | White | Micro Brush, 7 French diameter, 10" overall length | 3 | 23.00/pkg |
| B-8FR10 | Brown | Micro Brush, 8 French diameter, 10" overall length | 3 | 25.50/pkg |
| B-9FR10 | Red | Micro Brush, 9 French diameter, 10" overall length | 3 | 22.50/pkg |
| B-10FR10 | Green | Micro Brush, 10 French diameter, 10" overall length | 3 | 22.50/pkg |
| B-12FR10 | Black | Micro Brush, 12 French diameter, 10" overall length | 3 | 22.00/pkg |
| Optional Brushes: | | | | |
| B-8FR25 | Brown | Micro Brush, 8 French diameter, 25" overall length | 3 | \$25.50/pkg |
| B-11FR10 | Yellow | Micro Brush, 11 French diameter, 10" overall length | 3 | 22.50/pkg |
| B-14FR12 | Green | Micro Brush, 14 French diameter, 12" overall length | 3 | 22.50/pkg |
| B-14FR16 | Green | Micro Brush, 14 French diameter, 16" overall length | 3 | 23.00/pkg |
| B-14FR24 | Green | Micro Brush, 14 French diameter, 24" overall length | 3 | 23.50/pkg |
| B-20FR16 | White | Micro Brush, 20 French diameter, 16" overall length | 3 | 22.50/pkg |
| B-20FR24 | White | Micro Brush, 20 French diameter, 24" overall length | 3 | 24.00/pkg |
| B-29FR16 | Black | Micro Brush, 29 French diameter, 16" overall length | 3 | 22.50/pkg |
| B-29FR24 | Black | Micro Brush, 29 French diameter, 24" overall length | 3 | 23.50/pkg |
| B-2782 | n/a | French Scale Measuring Gauge FREE with first order of Micro Nest | 1 | 5.00 |

Minimum Order: \$25

Freight: F.O.B. Origin. Freight prepaid and added to invoice.

Shipping: UPS Ground Service, unless otherwise requested.

Terms: Net 30 days

Warranty: The ISI Micro Nest wall organizer is guaranteed for 90 days against manufacturer's defects. The brushes are NOT covered under the 90 day warranty.