

# ISI Instrument Care Solutions – Quick Reference Guide

Effective, biodegradable and non-hazardous • Safe for plastic, rubber and metal

ISI Product	Ready to use	Concentrate†	For ultrasonic cleaner	For washer sterilizer	Manual cleaning	Soak
<b>Enzymatic Detergents</b>						
Hemaway Enzymatic Original		●	●	●	●	●
Hemaway Enzymatic Tri-Max		●	●	●	●	●
Hemaway Foam	●				●	●
EndoTerg, for flexible scopes		●		●	●	
<b>Detergents</b>						
MPD, multi-purpose		●	●		●	●
Suds-Rite, low foam		●	●	●		
DuoKlenz, detergent, lubricant*		●	●	●		
<b>Rust &amp; Stain Removers</b>						
SSK-80, stainless steel cleaner		●	●	●		
CSK-81, cool soak		●			●	●
<b>Lubricants</b>						
Instra-Glide Soak	●					●
Instra-Glide Jet		●		●		
DuoKlenz, lubricant, detergent*		●	●	●		

\* DUAL PURPOSE CLEANER † CONCENTRATED FORMULA MAKES 128 GALLONS OF SOLUTION

## Free Samples

For more information,  
call 1-800-537-1945 or visit [www.isisurgery.com](http://www.isisurgery.com)

ISI Medical & O.R. Products [www.isisurgery.com](http://www.isisurgery.com) 1.800.537.1945