

Lubricants

Effective, biodegradable and non-hazardous
Safe for plastic, rubber, metal

Instra-Glide Soak

Non-silicone lubricant soak for surgical instruments

Instra-Glide Jet **NEW**

Water-soluble, non-silicone lubricant

Duo-Klenz **NEW**

Removes blood and organic materials while lubricating surgical instruments

Rust & Stain Removers

Effective, biodegradable and non-hazardous
Safe for plastic, rubber, metal

SSK-80

Stainless steel rust and stain remover for surgical instruments

CSK-81

Ready-to-use soak for pre-cleaning surgical instruments

Request Your
Free Samples

1-800-537-1945

www.isisurgery.com

Brushes, Sets & Organizers



ISI Brush Nest

The ISI Brush Nest is the wall organizer with a complete set of cleaning brushes in a variety of sizes and diameters.

Instrument Cleaning Brushes

ISI offers instrument cleaning brushes in a variety of lengths, diameters and compositions. Sold in sets and individually.



ISI Micro Nest

The ISI Micro Nest is the wall organizer with a complete set of micro brushes in a variety of sizes and diameters.

Micro Brushes

ISI Micro Brushes are available in a variety of French diameters. Sold in sets and individually.

 **Instrument Specialists, Inc.**

www.isisurgery.com

Toll Free in the USA: (800) 537-1945

Tel: (830) 249-9535 ▲ Fax: (830) 249-9433

32390 IH-10 West, Boerne, Texas 78006-9214

ISI Surgical Instrument Cleaners and Lubricants, Cleaning Brushes & Brush Organizers 2007



1-800-537-1945 ▲ www.isisurgery.com



**ISI Medical &
O.R. Products**

A Division of Instrument Specialists, Inc.

ISI Instrument Detergents, Cleaners, and Lubricants

Get the most from your surgical instruments with ISI instrument care solutions – FREE samples available



Back in the late 70's, our company's founder, Joseph R. Kabanek, Ph.D., noticed that many of the instruments sent for repair were not broken, but simply needed a thorough cleaning. His observation led to the development of our instrument care solutions – enzymatic detergents, detergents, rust and stain removers, and lubricants.

Today ISI continues to develop and perfect instrument care solutions to help keep delicate surgical and medical instruments functioning properly in hospitals across the country and around the world.

For your FREE samples, simply call our toll-free number or visit us at www.isisurgery.com

"...not broken, but simply needed a thorough cleaning..."



ISI Quick Reference Guide – Effective, biodegradable, non-hazardous and safe for plastic, rubber and metal

ISI Product	Ready to use	Concentrate†	For ultrasonic cleaner	For washer sterilizer	Manual cleaning	Soak
Enzymatic Detergents						
Hemaway Enzymatic Original		●	●	●	●	●
Hemaway Enzymatic Tri-Max		●	●	●	●	●
Hemaway Foam	●				●	●
EndoTerg, for flexible scopes		●		●	●	●
Detergents						
MPD, multi-purpose		●	●		●	●
Suds-Rite, low foam		●	●	●		
DuoKlenz, detergent, lubricant*		●	●	●		
Rust & Stain Removers						
SSK-80, stainless steel cleaner		●	●	●		
CSK-81, cool soak		●			●	●
Lubricants						
Instra-Glide Soak	●					●
Instra-Glide Jet		●		●		
DuoKlenz, lubricant, detergent*		●	●	●		

* Dual Purpose Cleaner † Concentrated formula makes 128 gallons of solution

Enzymatic Detergents

Effective, biodegradable and non-hazardous
Safe for plastic, rubber, metal

- Hemaway Original**
Enzymatic strength. Removes organic materials
- Hemaway Tri-Max**
Contains the three strongest enzymes available for maximum cleaning
- Hemaway Enzymatic Foam**
The ready-to-use, pre-cleaning foam
- EndoTerg NEW**
Low-foaming enzymatic detergent for flexible endoscopes. Bacteriostatic. Free rinsing

Detergents

Effective, biodegradable and non-hazardous
Safe for plastic, rubber, metal

- Duo-Klenz NEW**
Removes blood and organic materials while lubricating surgical instruments
- MPD**
The multi-purpose detergent
- Suds-Rite**
Low-foaming detergent concentrate for use in washer sterilizers and pressurized water units

Request Your Free Samples
1-800-537-1945
www.isisurgery.com

