

## ISI Surgical Table That Best Suits Your Needs – Quick Reference Guide

Standard Feature	Advanced Universal K	Universal K	Carter Hand Surgery Table	Self-Supporting Arm Table	Ultra-Light	Free-Standing Arm Table
Lexan Top		●		●		●
Carbonfiber Top	●				●	
Polycarbonate Top			●			
Weight	16 lbs	20 lbs	18 lbs	18 lbs	9 lbs	25 lbs
Telescoping Leg	●	●	●	optional	optional	●
Side Rail Mounting	●	●	●	●	●	●
Fluoroscopy	●	●	●	●	●	●
Traction	optional	optional	standard	optional	optional	optional
Under-Pad Mounting Piece	●	●	●	●	●	●
2" & 4" Foam Pads	optional	optional	optional*	optional	optional	optional
Height Adjustable Support Leg	●	●	●	optional	optional	●
Free Evaluation	●	●	●	●	●	●

\*The Carter Hand Surgery Table comes standard with a one-inch foam pad.

### Free Evaluations

For more information on ISI's 30 day, no-charge evaluations in your O.R., call 1-800-537-1945 or visit [www.isisurgery.com](http://www.isisurgery.com)

ISI Medical & O.R. Products [www.isisurgery.com](http://www.isisurgery.com) 1.800.537.1945