



ISI Medical & O.R. Products

A Division of Instrument Specialists, Inc. ▲ ISI for Ideas, Service, Integrity



Specialized Surgical Tables ▲ Alumi-Hand ▲ Digi-Grip
Finger Traps ▲ Brush Nest ▲ Micro-Brushes
Instrument Cleaning Solutions & Lubricants ▲ 1-800-537-1945

ISI Advanced Universal K Table

The Versatile Surgical Table for Hand Surgeries, Pediatric, Hand, Neck and Lower Extremity Surgeries



Combines Strength with Flexibility – and Weighs Only 16 lbs

The ISI Advanced Universal K Table is engineered to take advantage of the best materials available today: carbon-fiber composites, stainless steel and aluminum. The result is a portable surgical platform with unsurpassed durability weighing only 16 lbs!

The carbon-fiber composite table surface is radiolucent and the hourglass shape accommodates easy access for surgeons, staff, and equipment. Best yet, is the ability to use fluoroscopy without moving the patient; the entire table top allows fluoroscopy.

How Versatile?



The ISI Advanced Universal K Table is the ideal table for upper extremities, pediatric orthopaedic, hand and neck, and lower extremity surgeries. Plus, the optional traction package gives you the ability to perform traction cases.

Options Expand Upon Table's Versatility

▲ Traction Package



The Traction Package gives you the ability to perform traction surgeries on the upper extremity. The package includes a traction bar, four stainless steel finger traps with chain, four "S" hooks, 2 one-pound weights, 2 two-poundweights, and 1 five-pound weight in a durable carrying case.

▲ Under-Pad Mounting Piece



The Under-Pad Mounting Piece is used when side rails aren't available or when the Under-Pad Mounting Piece is preferred. Ideal when a stretcher is used.

The Under-Pad Mounting Piece replaces the standard mounting bracket of the ISI Advance Universal K Table. It simply slides under the stretcher's pad or O.R. table and is held in place by the patient's weight.

▲ 4" Foam Pad



The Advanced Universal K Table has a standard 2" foam pad to create a level surgical field. However, a 4" pad is an option when extra-thick O.R. table pads are used.

Weights Only 16 lbs



Free Evaluation

The ISI Advanced Universal K Table is available for 30 day, no-charge evaluations to allow you to try it in your O.R.

ISI Ultra-Light Specialty Table

The Legless Surgical Table Designed Specifically for Upper Extremity Surgeries



Portable and Light, this Surgical Table Stands Alone

Sturdy and “legless” – that’s the ISI Ultra-Light Specialty Table. Made of lightweight, carbon-fiber composites, this table is strong and yet weighs only 9 lbs. Every aspect of the ISI Ultra-Light Specialty Table has been designed and engineered to make it a top-of-the-line surgical platform.

ISI Ultra-Light Specialty Table offers the ultimate accessibility for surgeons, technicians and equipment during upper extremity surgeries. Convenient built-in stainless steel mounting clamps quickly attach to the O.R. table’s side rail. The carbon-fiber surgical surface allows for fluoroscopy without moving the patient.

For those times additional support is needed, ISI offers an optional leg. When not in use, the telescoping leg simply locks under the tabletop.

Weighs Only 9 lbs

Free Evaluation

The ISI Ultra-Light Specialty Table is available for 30 day, no-charge evaluations to allow you to try it in your O.R.

Options Expand Convenience Even Further

▲ Optional Leg



The telescoping Optional Leg may be used with the ISI Ultra-Light Specialty Table to provide even greater stability. Leg height is adjustable.

▲ 4" Foam Pad




The ISI Ultra-Light Specialty Table has a standard 2" foam pad to create a level surgical field. However, a 4" pad is available when extra-thick O.R. pads are used.

Compare ISI Specialty Tables

ISI Advanced Universal K Table for lower and upper extremity surgeries



ISI Ultra-Light Specialty Table for upper extremity surgeries



ISI Alumi-Hand

Surgical Hand Immobilizer – Safe, Sterile, Ready to Use



The ISI Alumi-Hand provides all the benefits of the Lead Hand and more.



With the Alumi-Hand you trade the hazards of lead for a sterile, disposable, affordable surgical hand immobilizer that allows for the use of X-ray and fluoroscopy! The Alumi-Hand is made with a pliable aluminum interior for easy hand positioning and a foam exterior for patient

comfort. Available in large and extra-large sizes.

Free Samples

Call ISI for your free Alumi-Hand samples.

1-800-537-1945



Customer Service

(800) 537-1945

ISI Digi Grip & Finger Traps

Upper Extremity Traction



The economical Digi Grip works with your existing traction apparatus and IV poles.

Digi Grip – Equipment with Lots of Pull!

The Digi Grip looks small but has all the strength you need for secure extremity traction – making it ideal for wrist arthroscopies and casting procedures. Constructed of powder-coated aluminum with 80-pound test stainless steel chain, the Digi Grip has 5 clips that make changing finger traps a snap. The economical Digi Grip works with your existing traction apparatus and IV poles. The ISI Digi Grip is proof that you can have both economy and quality.

ISI Stainless Steel and Nylon Finger Traps



These hard-to-find items are ideal for providing upper extremity traction and can be used with the ISI Digi Grip or ISI Advanced Universal K Table to create a complete traction system.


All traps come in a variety of sizes: extra-small, small, medium, large sizes, extra-large sizes. Available individually or in sets.

Stainless Steel Finger Traps are reusable and grip the patient's finger snugly. Wrap the patient's finger with sterile tape before application and the Stainless Steel Finger Trap will be extra secure without causing skin trauma.

Nylon Finger Traps are for optimum patient comfort. The Nylon Finger Trap are gentle on the patient's skin with a firm grip. Nylon Finger Traps are available in sterile and non-sterile packs.

ISI Brush Nest & Micro-Brushes

Infection Control Begins with Clean Instruments



The Brush Nest makes it easier for technicians to select the right brush for the cleaning task


ISI Brush Nest — Get Organized for Effective Instrument Cleaning

The ISI Brush Nest is a wall organizer with a complete set of cleaning brushes in a variety of sizes and diameters. Made of tough and durable plastic with a convenient perforated shelf, the Brush Nest is the ultimate organizer.

The Right Brush for the Right Task

The ISI Brush Nest comes complete with 18 specialized brushes in different lengths, diameters and composition. The long, thin brushes are perfect for cannulas and tubing. The smallest brushes are the answer for micros and the brass and stainless steel brushes make quick work of box locks. Individual packages of replacement brushes are available for quick delivery.

ISI Micro Brushes



These small brushes from ISI are designed specifically for those tiny, hard-to-clean instruments such as cannulas and are available in a variety of French diameters. Made of stain-

less steel handles and nylon bristles, the ISI Micro Brushes will last through countless cleanings.

ISI Instrument Cleaners & Lubricants

Effective, Biodegradable, and Non-Hazardous



Rust, Dried Blood, Organic Matter, Mineral Deposits and Stains.

These are the reasons why ISI offers a full line of Instrument Cleaners. ISI cleaning solutions give your technicians the means to help rid your instruments of debris, control cross-contamination and help maintain optimum instrument function. Here's the ISI lineup.

Enzymatic Detergents

Hemaway Enzymatic Foam – first step in cleaning
Hemaway Tri-Max – enzymatic detergent for blood & tissue removal

Detergents

MPD – multi-purpose detergent
Suds-Rite – low sudsing detergent

Rust & Stain Removers

SSK-80 – stainless steel cleaner for use directly in the ultrasonic
CSK-81 – cool soak solution for soaking instruments in basin

Lubricants

Insta-Glide – a non-silicone lubricant
Insta-Glide Target – an aerosol non-silicone “spot” lubricant

Call for Free Samples

Once you use ISI cleaning solutions and lubricants, we're sure you'll see the difference in your surgical instruments.

ISI for Ideas, Service, Integrity

Since 1978, the professionals at ISI have worked together with pride. A unique dedication to craftsmanship and quality at the best possible prices has earned ISI trust across the medical community.

Two Divisions, the Same Commitment

ISI Medical & O.R Products provides hand surgery accessories, specialty surgical tables, and a complete line of instrument cleaners and lubricants.

ISI Medical Equipment Service & Sales repairs and refurbishes rigid and flexible endoscopes, power and electronic equipment, and phaco handpieces. Because our technicians serve hospitals nationwide, ISI is an excellent source for refurbished equipment.



ISI Medical & O.R. Products

- ▲ Specialized Surgical Tables
- ▲ Alumi-Hand
- ▲ Dura-Hand
- ▲ Digi-Grip
- ▲ Finger Traps
- ▲ Brush Nest
- ▲ Micro Nest
- ▲ Micro-Brushes
- ▲ Instrument Cleaning Solutions & Lubricants



ISI Medical Equipment Service & Sales

Repairs for:

- ▲ Flexible & Rigid Endoscopes
- ▲ Pneumatic Handpieces
- ▲ Electrical Handpieces
- ▲ Surgical Cameras
- ▲ Phaco Handpieces
- ▲ Video Equipment
- ▲ Surgical Instruments
- ▲ Sales & Purchase of Refurbished Equipment

ISI, Always at Your Service

We look forward to the opportunity to serve our clients as well as earn the respect of new customers. For more information or to speak with one of our customer representatives, please call **1-800-537-1945**.



ISI Instrument Specialists, Inc.

Tel: (800) 537-1945 ▲ Fax: (830) 249-9433 ▲ www.isisurgery.com
32390 IH-10 West ▲ Boerne, Texas 78006-9214