

Tilt-Top Hand Table



ISI's Tilt-Top Hand Table is ideal for therapy, minor procedures, and examinations of the hand and upper extremity. Features include:

- Table top adjusts in height or angle, accommodating both right-handers or lefties
- Split-top allows larger top to be angled, while smaller top holds instruments, therapy items, or paperwork
- Integrated spring in leg allows for easy height adjustment
- Locking casters provide smooth mobility
- Cherry finish table and powder-coated base
- Removable pads for both sections of table top

Specifications:

Large work area: 21" x 17" x ¾"

Small work area: 9 ¾" x 16"

Base is 22 ¾" w x 17 ½" d and needs 3 ¼" clearance

Height adjustment: 26-37"

Call today for more information.

1.800.537.1945



ISI Medical & O.R. Products

Division of Instrument Specialists, Inc.

1.800.537.1945 FAX: 830.249.9433 www.isisurgery.com
32390 IH-10 West, Boerne, Texas 78006-9214