

# VBM Digital Tourniquets®

## Product Details

### Digital Display

- Large color coded LED display is easy to read
- Precise monitoring shows actual pressure cuff

### Precision Pressure Regulator

- Fast cuff inflation and deflation
- Regulator is easier to handle than touch pads
- Limits pressure to max. 600mmHg

### Flush Button

- To check for bleeding
- To release drug slowly after Bier's Block

### Alarm System

- Visual and audible alarm system to detect leaks, pressure deviations, and low battery voltage.

### Color Coded Extension Hoses

- With Positive Locking Connectors (PLC)
- Easy attachment and detachment
- For safe and leak-free inflation
- Kink resistant

### Automatic Timer

- Automatic timer recording after cuff inflation
- Provides elapsed cuff inflation time
- Audible signal every 30 minutes

### Battery Back-Up

- Long lasting Ni-MH Battery

### Built-in Data Port

- Supports optional printer for patient record
- Reports cuff pressure adjustments and total elapsed cuff inflation time on a self-adhesive label

### Stand with Basket

- Very stable
- Easy patient transport

**Call today for more information.**

**1.800.537.1945**







**ISI Medical & O.R. Products**




Division of Instrument Specialists, Inc.

1.800.537.1945 FAX: 830.249.9433 [www.isisurgery.com](http://www.isisurgery.com)  
32390 IH-10 West, Boerne, Texas 78006-9214

# Single Cuffs

Patient, Size	 <b>Fabric Cuffs</b> <ul style="list-style-type: none"> <li>Washable and reusable</li> <li>Fixation via Velcro tape</li> <li>Bladder and fabric sleeve can be replaced</li> </ul>			 <b>Silicone Cuffs</b> <ul style="list-style-type: none"> <li>Autoclavable at 134°C (273°F)</li> <li>Fixation via strap and buckle</li> <li>Easy to clean and dry quickly</li> </ul>			 <b>Silicone/Velcro Cuffs</b> <ul style="list-style-type: none"> <li>Autoclavable at 134°C (273°F)</li> <li>Fixation via Velcro tape</li> <li>Velcro sleeve can be replaced</li> <li>Easy to clean and dry quickly</li> </ul>			 <b>Disposable Cuffs</b> <ul style="list-style-type: none"> <li>Disposable, sterile peel pouch</li> <li>Low profile, welded cuff (no bladder)</li> <li>Eliminate risk of cross contamination</li> </ul>	
							PLC Connector	Luer Lock Connector			
Infant, 20cm (8")	20-54-700	—	20-64-700	20-34-700SLZ-1	20-34-700LL-1						
Child, 30cm (12")	20-54-710	20-94-710	20-64-710	20-34-710SLZ-1	20-34-710LL-1						
Arm, 35cm (14")	20-54-711	20-94-711	20-64-711	20-34-711SLZ-1	20-34-711LL-1						
Arm, 46cm (18")	20-54-712	20-94-712	20-64-712	20-34-712SLZ-1	20-34-712LL-1						

# Double Cuffs

Patient, Size	 <b>Fabric Cuffs</b> <ul style="list-style-type: none"> <li>Washable and reusable</li> <li>Fixation via Velcro tape</li> <li>Bladder and fabric sleeve can be replaced</li> </ul>			 <b>Silicone Cuffs</b> <ul style="list-style-type: none"> <li>Autoclavable at 134°C (273°F)</li> <li>Fixation via strap and buckle</li> <li>Easy to clean and dry quickly</li> </ul>			 <b>Silicone/Velcro Cuffs</b> <ul style="list-style-type: none"> <li>Autoclavable at 134°C (273°F)</li> <li>Fixation via Velcro tape</li> <li>Velcro sleeve can be replaced</li> <li>Easy to clean and dry quickly</li> </ul>		
Infant, 20cm (8")	20-50-700	—	20-60-700						
Child, 30cm (12")	20-50-710	20-90-710	20-60-710						
Arm, 35cm (14")	20-50-711	20-90-711	20-60-711						
Arm, 46cm (18")	20-50-712	20-90-712	20-60-712						

Call today for more information.  
**1.800.537.1945**



**ISI Medical & O.R. Products**

Division of Instrument Specialists, Inc.

1.800.537.1945 FAX: 830.249.9433 [www.isisurgery.com](http://www.isisurgery.com)  
 32390 IH-10 West, Boerne, Texas 78006-9214