

VBM Roll-On Cuff



- Allows fast and easy blood exsanguinations within seconds (as alternative to Esmarch bandage)
- Ideal for hand surgery
- Material: Silicone, transparent
- Reusable, autoclavable up to 134°C (273°F)
- Latex-free



Exsanguinator

Fast and easy blood exsanguinations prior to surgery with single cuffs. Easily and smoothly rolls over the limb even with I.V. Catheter in place (for Bier's Block).



Exsanguinator and Tourniquet

If used to maintain a bloodless field. The brake block keeps the cuff in place and prevents it from rolling back.

Call today for
more information.
1.800.537.1945



ISI Medical & O.R. Products

Division of Instrument Specialists, Inc.

1.800.537.1945 FAX: 830.249.9433 www.isisurgery.com
32390 IH-10 West, Boerne, Texas 78006-9214

Roll-On Cuffs

- Exsanguinator
allows fast and easy blood exsanguination within seconds (as alternative to Esmarch Bandage)
- Material: Silicone, transparent
- Reusable, autoclavable up to 134°C (273°F)
- Latex free

Size	Color code	Patient limb		REF
0	grey	16-21cm	6"-8"	21-95-710
1	white	20-28cm	8"-11"	21-95-711
2	light-blue	26-33cm	10"-13"	21-95-712
3	yellow	31-40cm	12"-16"	21-95-719
4	red	38-49cm	15"-19"	21-95-722
5	green	47-60cm	18"-23"	21-95-727
6	dark-blue	58-70cm	23"-28"	21-95-729



Handpump

with manometer and
1m (40") connecting tube

REF 21-95-350



Manometer

with 1m (40") connecting tube

REF 21-95-300

Brake-Block

made of Silicone,
autoclavable up to 134°C (273°F)

REF 21-95-200 small size
REF 21-95-201 large size



Tape-Measure

color coded, for easy identification
of patient extremity and cuff size

REF 21-95-111



Roll-On Cuff Set – small sizes

complete case with Roll-On Cuff sizes 1-4, manometer with
connecting tube, tape-measure, brake-block small and large

REF 21-95-400

Roll-On Cuff Set – large sizes

complete case with Roll-On Cuff sizes 3-6, handpump with
manometer and connecting tube, tape-measure, brake-block
small and large

REF 21-95-401

Call today for
more information.
1.800.537.1945



ISI Medical & O.R. Products

Division of Instrument Specialists, Inc.

1.800.537.1945 FAX: 830.249.9433 www.isisurgery.com
32390 IH-10 West, Boerne, Texas 78006-9214