


ISI Surgical Instrument Cleaners and Lubricants, Cleaning Brushes and Brush Organizers



A Division of Instrument Specialists, Inc.

2007

1-800-537-1945 www.isisurgery.com
ISI for Ideas, Service, Integrity

Enzymatic Detergents
Detergents
Rust & Stain Removers
Lubricants
Instrument Brushes
Micro Brushes
ISI Brush Set Organizers

ISI Instrument Detergents, Cleaners, and Lubricants

Get the most from your surgical instruments with ISI instrument care solutions – FREE samples available



Back in the late 70's, our company's founder, Joseph R. Kabanek, Ph.D., noticed that many of the instruments sent for repair were not broken, but simply needed a thorough cleaning. His observation led to the development of our instrument care solutions – enzymatic detergents, detergents, rust and stain removers, and lubricants.

Today ISI continues to develop and perfect instrument care solutions to help keep delicate surgical and medical instruments functioning properly in hospitals across the country and around the world.

For your FREE samples, simply call our toll-free number or visit us at www.isisurgery.com

"...not broken, but simply needed a thorough cleaning..."



Free Samples 1-800-537-1945 www.isisurgery.com

ISI Quick Reference Guide

Effective, biodegradable and non-hazardous • Safe for plastic, rubber and metal

ISI Product	Ready to use	Concentrate [†]	For ultrasonic cleaner	For washer sterilizer	Manual cleaning	Soak
Enzymatic Detergents						
Hemaway Enzymatic Original		●	●	●	●	●
Hemaway Enzymatic Tri-Max		●	●	●	●	●
Hemaway Foam	●				●	●
EndoTerg, for flexible scopes		●		●	●	●
Detergents						
MPD, multi-purpose		●	●		●	●
Suds-Rite, low foam		●	●	●		
DuoKlenz, detergent, lubricant*		●	●	●		
Rust & Stain Removers						
SSK-80, stainless steel cleaner		●	●	●		
CSK-81, cool soak		●			●	●
Lubricants						
Instra-Glide Soak	●					●
Instra-Glide Jet		●		●		
DuoKlenz, lubricant, detergent*		●	●	●		

* Dual Purpose Cleaner † Concentrated formula makes 128 gallons of solution

■ Enzymatic Detergents

Effective, biodegradable and non-hazardous

Hemaway Original

Enzymatic strength removes organic materials

Hemaway Original removes dried blood, tissue, mucous, and other body fluids. Hard to reach places like box joints and tubing can't hide from the cleaning power of Hemaway Original. Reduces scrubbing time and prolongs the life of your instruments, Hemaway Original is formulated for daily use.

- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Removes organic materials, even dried-on blood
- ▶ Safe for plastic, rubber, metal

Hemaway Tri-Max

Hemaway Tri-Max contains the three strongest enzymes available for maximum cleaning

Similar to Hemaway Original, but formulated with the strongest combination of enzymes available. Concentrated Hemaway Tri-Max greatly reduces the need for scrubbing and prolongs the life of your instruments. Hemaway Tri-Max is strong enough to remove dried blood, tissue, mucous, and other body fluids, yet gentle enough to be used everyday.

- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Removes organic materials, even dried-on blood
- ▶ Maximum combination of enzymes available
- ▶ Safe for plastic, rubber, metal

Hemaway Enzymatic Foam

The ready-to-use, pre-cleaning foam for surgical instruments

Hemaway Enzymatic Foam is ideal for use in the O.R. to pre-clean surgical instruments and prevent corrosion. The foam contains active enzymes that go to work immediately to break down organic materials.

- ▶ Ready to use; spray directly on surgical instruments
- ▶ Removes organic materials safely
- ▶ Safe for plastic, rubber, metal

EndoTerg **NEW**

The low-foaming enzymatic detergent for flexible endoscopes

EndoTerg breaks down dried blood, tissue, mucous, and other bio-burdens on flexible endoscopes without leaving a soapy residue. Designed specifically for flexible endoscopes, EndoTerg provides thorough but gentle cleaning.

- ▶ Bacteriostatic and free rinsing
- ▶ Removes organic materials safely
- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Safe for plastic, rubber, metal

■ Detergents

Effective, biodegradable and non-hazardous

Duo-Klenz **NEW**

Removes blood and organic materials while lubricating surgical instruments

Highly concentrated, low-foaming Duo-Klenz removes blood and organic materials while safely lubricating surgical instruments – all in one easy step! Duo-Klenz is formulated for use in the wash cycle of washer sterilizers, cart washers, bottle washers, and other pressurized water units.

- ▶ Detergent combined with a non-silicone lubricant
- ▶ Removes organic materials, even dried-on blood
- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Safe for plastic, rubber, metal

MPD

The multi-purpose detergent

MPD is the highly concentrated multi-purpose detergent formulated to clean surgical instruments manually or with an ultrasonic cleaner. MPD is tough enough to remove blood and organic material and safe enough for everyday cleaning.

- ▶ For manual cleaning and use in ultrasonics
- ▶ Removes organic materials, even dried-on blood
- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Safe for plastic, rubber, metal

Suds-Rite

Low-foaming detergent concentrate for use in washer sterilizers and other pressurized water units

Use Suds-Rite in the wash cycle of washer sterilizers, cart washers, bottle washers, and other pressurized water units. Suds-Rite removes blood and organic materials from surgical instruments.

- ▶ The anti-foaming agent makes it perfect for use in the wash cycle of washer sterilizers and pressurized water units
- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Removes organic materials, even dried-on blood
- ▶ Safe for plastic, rubber, and metal



■ Rust & Stain Removers

Effective, biodegradable and non-hazardous

SSK-80

The stainless steel rust and stain remover for surgical instruments

SSK-80 is our best and most convenient instrument rust and stain remover for use in ultrasonic cleaners.

Used routinely, SSK-80 will extend the life of your instruments by removing and inhibiting rust, stains, and mineral deposits. The result: pivot points, box joints, and moving instrument parts function better. Simply add SSK-80 to your ultrasonic during the fill cycle. It's that easy!

- ▶ Concentrated; pour directly into the ultrasonic cleaner
- ▶ Removes rust, stains and mineral deposits
- ▶ Does not leave a residue
- ▶ Safe for plastic, rubber, metal

CSK-81

Concentrated rust and stain remover for soaking surgical instruments

Don't have an ultrasonic cleaner available? Or need to clean instruments not suitable for ultrasonic cleaners? Either way, CSK-81 instrument soaking solution removes rust, stains and mineral deposits leaving your surgical instruments cleaner and in better condition. Use routinely to inhibit new rust and extend the life of your instruments. Simply soak your surgical instruments and CSK-81 will do the rest!

- ▶ Concentrated formula makes 128 gallons of solution
- ▶ No heating required
- ▶ Removes rust, stains and mineral deposits
- ▶ Safe for plastic, rubber, metal



■ Lubricants

Effective, biodegradable and non-hazardous

Instra-Glide Soak

Non-silicone lubricant soak for surgical instruments

Instra-Glide Soak is the non-silicone lubricant that protects your surgical instruments. Instra-Glide Soak is water-soluble and provides a protective film without getting sticky or gummy. Keep your surgical instruments functioning at their optimum with this non-silicone lubricant.

- ▶ Ready to use; soak instruments prior to sterilization
- ▶ Water-soluble for easy rinsing
- ▶ Safe for plastic, rubber, metal

Instra-Glide Jet **NEW**

Non-silicone lubricant

Instra-Glide Jet is the non-silicone lubricant designed to keep surgical instruments functioning properly. Water-soluble Instra-Glide Jet provides a protective film without getting sticky or gummy. Helps prevent rust, stains and corrosion.

- ▶ Concentrated – makes 128 gallons of solution
- ▶ Water-soluble for easy rinsing
- ▶ Preservatives added to prolong shelf life
- ▶ Use in last cycle of washer to lubricate set

Duo-Klenz **NEW**

Removes blood and organic materials while lubricating surgical instruments

Highly concentrated, low-foaming Duo-Klenz removes blood and organic materials while safely lubricating surgical instruments – all in one easy step! Duo-Klenz is formulated for use in the wash cycle of washer sterilizers, cart washers, bottle washers, and other pressurized water units.

- ▶ Detergent combined with a non-silicone lubricant
- ▶ Removes organic materials, even dried-on blood
- ▶ Concentrated formula makes 128 gallons of solution
- ▶ Safe for plastic, rubber, metal



**Request Your
Free Samples**
1-800-537-1945
www.isisurgery.com

Brush Sets & Organizers

ISI Brush Nest, ISI Micro Nest wall-mount organizers

Get organized for effective instrument cleaning

The ISI Brush Nest and Micro Nest wall organizers come with a complete set of cleaning brushes in a variety of sizes and diameters. Both the Brush Nest and Micro Nest are designed to help technicians select the best brush for the cleaning task.

Individual packages of replacement brushes are available for both the ISI Brush Nest and ISI Micro Nest.



ISI Brush Nest

Made of tough durable plastic with a convenient self-draining shelf, the ISI Brush Nest comes with 18 specialized brushes in different lengths, diameters and compositions.



ISI Micro Nest

Made of durable plastic, the ISI Micro Nest organizes a complete set of micro brushes. All micro brushes have stainless steel handles and nylon bristles, plus come in a variety of sizes and French diameters. These brushes are designed for hard-to-clean instruments such as cannulas.



Brushes

ISI brushes are durable and come in a variety of diameters and sizes



Instrument Cleaning Brushes

ISI offers durable instrument cleaning brushes in a variety of lengths, diameters and compositions. Sold in sets and individually.

Micro Brushes

ISI Micro Brushes are available in a variety of French diameters. Made with stainless steel handles and nylon bristles, ISI Micro Brushes will last through countless cleanings. Sold in sets and individually.

ISI Free Offers

Call ISI at 1-800-537-1945 or visit www.isisurgery.com

Free samples of ISI instrument cleaners and lubricants

If you haven't tried our instrument cleaners and lubricants, here's your chance! Call or go online for your free samples.

Free Alumi-Hand samples

With the Alumi-Hand you trade the hazards of lead for a sterile, disposable, affordable surgical hand immobilizer. Call or go online for your free sample.

Free ISI surgical table evaluations

Don't just take our word that these are great surgical platforms – see for yourself! The Universal K, Advanced Universal K, Self-Supporting Arm Table, Ultra-Light, and Free-Standing Arm Table are available for a 30 day, no-charge evaluation in your O.R.

Two divisions, the same commitment

Since 1978, the professionals at ISI have worked together with pride. A unique dedication to craftsmanship and quality at the best possible price has earned ISI trust across the medical community.

ISI Medical & O.R. Products

- ▲ Specialty Surgical Tables
- ▲ Leadless Hand Immobilizers
- ▲ Upper Extremity Traction Devices
- ▲ Instrument Cleaning Solutions
- ▲ Instrument Lubricants
- ▲ Instrument Cleaning Brushes, Sets and Brush Organizers
- ▲ GSA's Federal Supply Service (FSS)
- ▲ TXMAS Contract Pricing
- ▲ Other GPO Pricing Available

ISI Medical Equipment Service & Sales

Service and Repairs:

- ▲ Flexible & Rigid Endoscopes
- ▲ Medical Camera Repair
- ▲ Monitor & Printer Repair
- ▲ Powered Hand Piece Repair
- ▲ Phaco Systems Repair
- ▲ Endoscopic Accessories
- ▲ GSA's Federal Supply Service (FSS)
- ▲ TXMAS Contract Pricing
- ▲ Other GPO Pricing Available

ISI, always at your service


We look forward to the opportunity to serve our clients as well as earn the respect of new customers. For more information visit www.isisurgery.com or speak with one of our customer representatives at 1-800-537-1945.

Instrument Specialists, Inc.

www.isisurgery.com ▲ Toll Free in the USA: (800) 537-1945

Tel: (830) 249-9535 ▲ Fax: (830) 249-9433

32390 IH-10 West, Boerne, Texas 78006-9214





O.R. & Medical Products

- ▲ Advanced Universal K Table
- ▲ ISI Ultra-Light Specialty Table
- ▲ Alumi-Hand
- ▲ Digi-Grip & Finger Traps
- ▲ Specialized Instrument Cleaning Solutions
- ▲ ISI Brush Nest
- ▲ Specialty Brushes

Medical Equipment Sales & Service

Repairs for:

- ▲ Flexible & Rigid Endoscopes
- ▲ Pneumatic Handpieces
- ▲ Electrical handpieces
- ▲ Surgical Cameras
- ▲ Phaco Handpieces
- ▲ Video Equipment
- ▲ Surgical Instruments
- ▲ Sales & Purchase of Refurbished Equipment

Detergents		
	Description	Price
MPD: Multi-Purpose Detergent		
M5100	1 gallon	\$18
M5110	case of 4 gallons	\$63
M5120	5 gallon keg	\$77
M5125	15 gallon keg	\$225
M5127	30 gallon keg	\$420
M5130	55 gallon barrel	\$545
M5140	Pallet of 48 cases	\$2,256
Suds-Rite: Low Sudsing Detergent for Pressurized Cleaning Units		
SR6100	case of 4 gallons	\$60
SR6120	5 gallon	\$65
SR6125	15 gallon	\$195
SR6130	30 gallon drum	\$353
SR6135	55 gallon barrel	\$495
Duo-Klenz: Combination low-sudsing detergent with non-silicone lubricant		
DK3100	case of 4 gallons	\$68
DK3120	5 gallon keg	\$74
DK3125	15 gallon keg	\$220
DK3130	30 gallon keg	\$399
DK3135	55 gallon barrel	\$560

Lubricants		
	Description	Price
Instra-Glide Soak: Instrument Lubricant		
IG8705	1 gallon	\$18
IG8715	case of 4 gallons	\$66
IG8720	5 gallon keg	\$78
IG8725	15 gallon keg	\$179
IG8727	30 gallon keg	\$351
Instra-Glide Jet: Instrument Lubricant for pressurized Cleaning Units		
J4100	1 gallon	\$17
J4110	case of 4 gallons	\$64
J4120	5 gallon keg	\$76
J4125	15 gallon keg	\$175
J4130	30 gallon keg	\$340
J4155	55 gallon barrel	\$552

Rust & Stain Removers		
	Description	Price
SSK-80: Rust & Stain Remover for use in the Ultrasonic		
S5200	1 gallon	\$21
S5210	case of 4 gallons	\$83
S5220	5 gallon keg	\$99
S5225	15 gallon keg	\$297
S5227	30 gallon keg	\$590
S5230	55 gallon barrel	\$786
S5240	Pallet of 48 cases	\$2,976
CSK-81: Rust & Stain Remover for Soaking Instruments		
C5300	1 gallon	\$23
C5310	case of 4 gallons	\$86
C5320	5 gallon keg	\$105
C5325	15 gallon keg	\$305
C5327	30 gallon keg	\$610
C5330	55 gallon barrel	\$847
C5340	Pallet of 48 cases	\$3,096


Enzymatic Detergents		
	Description	Price
Hemaway Original: Enzymatic Detergent		
H7005	1 gallon	\$18
H7010	case of 4 gal.	\$63
H7020	5 gallon keg	\$77
H7025	15 gallon keg	\$225
H7027	30 gallon keg	\$420
Hemaway Tri-Max: Enzymatic Detergent		
H7605	1 gallon	\$21
H7615	case of 4 gal.	\$81
H7620	5 gallon keg	\$86
H7625	15 gallon keg	\$239
H7627	30 gallon keg	\$404
Hemaway Enzymatic Foam:		
H7700	16 oz bottles/case of 6	\$52
EndoTerg Enzymatic Scope Detergent:		
ET-9500	1 gal.	\$21
ET-9515	case of 4 gal.	\$81
ET-9520	5 gal. Keg	\$86
ET-9525	15 gal. Keg	\$239
ET-9527	30 gal. Keg	\$404

ISI Brush Nest – Brush set and organizer for instrument cleaning brushes

The ISI Brush Nest is the wall organizer with a complete set of 18 instrument cleaning brushes in a variety of sizes, diameters and composition. Made of tough, durable plastic with a convenient perforated shelf, the ISI Brush Nest is the ultimate organizer. Technicians can get to work right away with the right brushes for the task.

The long, thin brushes are perfect for cannulas and tubing. The smallest brushes are the answer for micros. The brass and stainless steel brushes make quick work of box locks.

Individual packages of replacement brushes are available for quick delivery.



ISI Medical & O.R. Products

A DIVISION OF INSTRUMENT SPECIALISTS, INC.

1-800-537-1945 www.isisurgery.com

Stock No.	Color Code	Description of Overall Length and Diameter	# Per Pkg	Price Per Pkg
BN-9600		Brush Nest set – the wall hanger with 18 different brushes	1	\$235.00
BN-9650		Re-order set of one of each standard brush	18	62.00
Standard Brushes				
BN-130-1	White	21" tube brush, 20 French, 1/4" diameter with 4" bristle length	3	\$18.50/pkg
BN-130-3	Black	21" tube brush, 30 French, 3/8" diameter with 4" bristle length	3	18.50/pkg
BN-56740	Blue	20" tube brush, 24 French, 5/16" diameter with 3-1/2" bristle length	3	17.25/pkg
BN-660-1	Red	6-1/4" tube brush, 9 French, 1/8" diameter with 2-1/4" bristle length	6	16.25/pkg
BN-660-5	Blue	6-1/4" tube brush, 24 French, 5/16" diameter with 2-1/4" bristle length	6	16.25/pkg
BN-661	Yellow	5" tube brush, 21 French, 9/32" diameter with 1-1/2" bristle length	6	15.75/pkg
BN-662	Red	5" tube brush, 9 French, 1/8" diameter with 1-1/2" bristle length	6	15.75/pkg
BN-663	Green	5" tube brush, 14 French, 3/16" diameter with 1-1/2" bristle length	6	15.75/pkg
BN-670	White	6-1/4" tube brush, 20 French, 1/4" diameter with 2" bristle length	6	16.25/pkg
BN-671	Black	6-1/4" tube brush, 30 French, 3/8" diameter with 2" bristle length	6	16.25/pkg
BN-940-1	White	12" tube brush, 20 French, 1/4" diameter with 4" bristle length	3	13.85/pkg
BN-940-3	Black	12" tube brush, 30 French, 3/8" diameter with 4" bristle length	3	13.85/pkg
BN-SX3198	Gray	16" tube brush, 34 French, 1/2" diameter with 3" bristle length	3	15.25/pkg
BN-SX3923	Brown	13" tube brush, 34 French, 1/2" diameter with 3" bristle length	3	14.85/pkg
BN-069		Nylon brush – tooth brush style	12	25.00/pkg
BN-83112		Brass brush – tooth brush style	12	25.00/pkg
BN-83111		Stainless steel brush – tooth brush style	12	25.00/pkg
BN-3000		Nylon nail scrub brush	3	9.50/pkg
Optional Brushes				
BN-4100		Pipe cleaners; 6" long, 3/16" diameter	36	\$ 3.00/pkg
BN-4200		Pipe cleaners; 12" long, 3/16" diameter	24	4.00/pkg
BN-4300		Pipe cleaners; 12" long, 3/8" diameter	24	7.50/pkg
B-06213		Small trach brushes 3-3/4"; fits tube sizes 000-5; 1-7/8" long, 1/2" diameter	3	12.00/pkg
B-06214		Large trach brushes 8"; fits tube sizes 6-12; 3-3/4" long, 1/2" diameter	3	14.00/pkg
BN-130-18C	Red	21" tube brush, 9 French, 1/8" diameter with 2" bristle length	3	21.50/pkg
BN-56737	Blue	14" tube brush, 3/4" diameter with 3-1/2" bristle length	3	15.25/pkg
BN-56738	Yellow	14" tube brush, 7/8" diameter with 4" bristle length	3	15.25/pkg
BN-56739	Blue	8" tube brush, 3/4" diameter with 3-1/3" bristle length	6	19.85/pkg
BN-660-316	Green	6-1/4" tube brush, 14 French, 3/16" diameter with 2-1/4" brush length	6	16.25/pkg
BN-660-14	White	6-1/4" tube brush, 20 French, 1/4" diameter with 2-1/4" brush length	6	16.25/pkg

Minimum Order: \$25

Freight: F.O.B. Origin. Freight prepaid and added to invoice.

Shipping: UPS Ground Service, unless otherwise requested.

Terms: Net 30 days

Warranty: The ISI Brush Nest wall organizer is guaranteed for 90 days against manufacturer's defects. The brushes are not covered under the 90 day warranty.


ISI Micro Nest – Brush set and organizer for micro brushes

The ISI Micro Nest is the wall-mounted organizer with a complete set of micro brushes in a variety of sizes and diameters. The ISI Micro Nest makes it easier for technicians to select the right brush for the cleaning task.

ISI Micro Brushes are designed specifically for those hard-to-clean instruments such as

cannulas. ISI Micro Brushes are available in a variety of French diameters. Made with stainless steel handles and nylon bristles, the ISI Micro Brushes will last through countless cleanings.

Brush packages may be ordered separately.



ISI Medical & O.R. Products

A DIVISION OF INSTRUMENT SPECIALISTS, INC.

1-800-537-1945 www.isisurgery.com

Stock No.	Color Code	Description of Overall Length and Diameter	# Per Pkg	Price Per Pkg
M-2004	n/a	ISI Micro Nest set – the wall organizer with 8 micro brushes (one each listed below)	1	\$85.00
M-2800	n/a	Micro Nest Wall Organizer – the organizer without brushes	1	25.00
M-2650	n/a	Re-order set of one of each micro brush	9	62.00
Micro Nest Brushes:				
B-4FR10	Gray	Micro Brush, 4 French diameter, 10" overall length	3	\$23.50/pkg
B-5FR10	Yellow	Micro Brush, 5 French diameter, 10" overall length	3	23.50/pkg
B-6FR10	Blue	Micro Brush, 6 French diameter, 10" overall length	3	23.00/pkg
B-7FR10	White	Micro Brush, 7 French diameter, 10" overall length	3	23.00/pkg
B-8FR10	Brown	Micro Brush, 8 French diameter, 10" overall length	3	25.50/pkg
B-9FR10	Red	Micro Brush, 9 French diameter, 10" overall length	3	22.50/pkg
B-10FR10	Green	Micro Brush, 10 French diameter, 10" overall length	3	22.50/pkg
B-12FR10	Black	Micro Brush, 12 French diameter, 10" overall length	3	22.00/pkg
Optional Brushes:				
B-8FR25	Brown	Micro Brush, 8 French diameter, 25" overall length	3	\$25.50/pkg
B-11FR10	Yellow	Micro Brush, 11 French diameter, 10" overall length	3	22.50/pkg
B-14FR12	Green	Micro Brush, 14 French diameter, 12" overall length	3	22.50/pkg
B-14FR16	Green	Micro Brush, 14 French diameter, 16" overall length	3	23.00/pkg
B-14FR24	Green	Micro Brush, 14 French diameter, 24" overall length	3	23.50/pkg
B-20FR16	White	Micro Brush, 20 French diameter, 16" overall length	3	22.50/pkg
B-20FR24	White	Micro Brush, 20 French diameter, 24" overall length	3	24.00/pkg
B-29FR16	Black	Micro Brush, 29 French diameter, 16" overall length	3	22.50/pkg
B-29FR24	Black	Micro Brush, 29 French diameter, 24" overall length	3	23.50/pkg
B-27B2	n/a	French Scale Measuring Gauge FREE with first order of Micro Nest	1	5.00

Minimum Order: \$25

Freight: F.O.B. Origin. Freight prepaid and added to invoice.

Shipping: UPS Ground Service, unless otherwise requested.

Terms: Net 30 days

Warranty: The ISI Micro Nest wall organizer is guaranteed for 90 days against manufacturer's defects. The brushes are NOT covered under the 90 day warranty.